Tizian Blank NEM

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Instructions

Description / indications

Tizian Blank NEM are made of a biocompatible, millable non-precious alloy on CoCr-basis. Nickel- and beryllium-free.

Tizian Blank NEM are indicated for the production of crowns and bridges. Adherence oxid creator especially suited to the crown and bridge technique provide for an excellent metal ceramic bond. It is up to the attending dentist's decision whether the material should be used for special restoration designs such as free-end bridges or attachments.

Chemical composition

Co [% max.]	Cr [% max.]	W [% max.]	Mn [% max.]	Fe [% max.]	Si [% max.]	C [% max.]
remainder	28.0	8.5	0.25	< 0.5	1.65	< 0.1

General information

Type	4, dental alloy on CoCr-basis
Color	silver-colored
Joining technique	laser welding, soldering
Suitable ceramics	conventional ceramics, e. g. Nuance 850

Physical properties*

density	[g/cm ³]	8.3
E-module	[GPa]	190
Vicker's hardness	[HV10]	< 320
melting point (solidus/liquidus)	[°C]	1390 / 1415
coefficient of thermal expansion	[*10 ⁻⁶ K ⁻¹]	14.5
elongation at break	[%]	10

* The physical properties are average values and may deviate slightly from one lot to another.

Safety note:

When processing non-precious alloy blanks, dust is emitted which can damage the lungs, and irritate skin and eyes. Therefore, only process the material while running an adequate extractor system.

Please adhere to the following safety rules:

- Process non-precious alloy blanks only underneath a local extractor system.
- Ensure that product dust will not come in touch with eyes or mucous membranes.
- Do not inhale dust.
- Wear face mask and protective goggles.

1. Shaping (form grinding)

To manufacture and process dental restorations made of Tizian NEM non-precious alloy blanks, only use machines, tools and milling systems which are suited to handle non-precious alloy blanks. Please adhere to the instructions included with the machine, tool or milling system.

2. Framework dimensions

Dimension the frameworks as usual. Ensure stable interdental connections. Apply metal collars to the bridge units to avoid heat accumulation when firing a ceramic facing. For composite facings, prepare retentions on the framework.

3. Detaching frameworks from the blank

Detach the milled frameworks from the blank with help of hard metal mills or separating discs (0.3 mm recommended) that are suited to machining non-precious alloy blanks. Pay particular attention to the margins of the framework.

4. Finishing the frameworks

Use only clean hard metal mills for finishing, without applying pressure and at a maximum of 25,000 rpm.

Residual oxides can be removed with polishing beads. Plane the finished metal surface with ceramic bonded grinding tools.

For the polish, we recommend standard rubber polishers, polishing paste, brushes and cotton buffs (by Schütz Dental).

Clean frameworks with a steam cleaner or in an ultrasonic bath with distilled water. Do not immerse non-precious alloy objects in an etching bath.

5. Ceramic facing

Tizian NEM objects can be faced with all conventional high-fusing ceramics (e.g. Nuance 850) which exhibit a suitable thermal expansion range. The material is suited to temperatures up to 980 °C. Avoid exceeding this temperature.

An oxide firing is not compulsory, but it is advantageous. For glass ceramics, we recommend an oxide firing temperature of 900 °C. For classic feldspath ceramics with higher individual firing temperatures, an oxide firing temperature of up to 950 °C is acceptable. Generally, fire under vacuum for 10 min.

After the oxide firing, remove the oxide by blasting with aluminium oxide (max. 110 µm) at 2 bar and in a 45° angle.

Using a bonding agent for non-precious alloys will positively influence the metal-ceramic bond.

Depending on the size of the framework, we recommend cooling it down slowly.

Adhere to the instructions of the ceramic's manufacturer.

6. Joining technique

Objects made of Tizian NEM can be soldered or laser welded. We also recommend to forego the use of a post-solder. For laser welding, we recommend the use of e.g. CoCr laser welding wire by Schütz Dental.

7. Fixing techniques

Fix the restoration with conventional phosphate or glass ionomer cements or with an adhesive technique.

Only use materials suitable for use with non-precious alloys. Adhere to the manufacturer's instructions.

Storage:

Store blanks in a dry place. Generally avoid contamination.

Disposal:

Dispose of metal scraps in accordance with local pollution control regulations.

Important:

Do not reuse Tizian NEM.

Side-effects:

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (e.g. allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

Contra-indications / interactions:

If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. In such cases, we will supply the composition of our medical device upon request. The dentist should consider known interactions and cross-reactions of the product with other materials already in the patient's mouth before using the product.

Note:

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

Guarantees:

Our processing instructions, regardless of whether they are issued verbally, in writing or in the form of a practical demonstration, are based upon our own experience and can therefore be regarded solely as guidelines.

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Schütz Dental GmbH · Dieselstr. 5-6 · 61191 Rosbach/Germany