



SCHÜTZ-DENTAL GmbH

Zusammensetzungen der Vor- und Nachlote

Unsere Lote haben folgende Zusammensetzungen (in %):

Element	Au	Pd	Pt	Ag	Cu	Zn	In	Sonstige
Vorlot für								
Alphador Nr. 1	64,0	---	0,4	34,9	---	---	0,6	Ir 0,1
Alphador Inlay	68,0	---	---	13,3	12,8	4,9	---	Sn 1,0
Alphador LF	68,0	---	---	13,3	12,8	4,9	---	Sn 1,0
Alphador HP	68,0	---	---	13,3	12,8	4,9	---	Sn 1,0
Alphador PDF	60,0	---	0,5	36,4	---	2,0	---	Sn 1,0
Alphador Soleil	68,0	---	---	13,3	12,8	4,9	---	Sn 1,0
Alphador 515	86,2	5,0	5,0	---	---	3,0	0,8	
Alphador dent U	78,0	10,0	---	---	0,5	1,5	---	Ni 10,0
Elceral	86,2	5,0	5,0	---	---	3,0	0,8	
Elceral CF	86,2	5,0	5,0	---	---	3,0	0,8	
Nachlot für								
Alphador Inlay	59,9	---	0,2	14,6	15,2	4,0	6,0	Ir 0,1
Alphador LF	59,9	---	0,2	14,6	15,2	4,0	6,0	Ir 0,1
Alphador HP	59,9	---	0,2	14,6	15,2	4,0	6,0	Ir 0,1
Alphador PDF	42,0	---	---	41,3	2,0	11,0	3,7	
Alphador Soleil	59,9	---	0,2	14,6	15,2	4,0	6,0	Ir 0,1
Alphador 515	42,0	---	---	41,3	2,0	11,0	3,7	
Alphador dent U	42,0	---	---	41,3	2,0	11,0	3,7	
Elceral	42,0	---	---	41,3	2,0	11,0	3,7	
Elceral CF	42,0	---	---	41,3	2,0	11,0	3,7	

Die Mengenangaben der Bestandteile sind Mittelwerte und können von der aktuellen Liefercharge geringfügig abweichen - nach DIN EN 29333 max. 0,5 %.

Rosbach, den 17.04.2000